

Zinc finger matrin-type protein 3

Catalog No: #AP89206



Package Size: #AP89206-1 50ug #AP89206-2 100ug #AP89206-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Zinc finger matrin-type protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Wild-type p53-induced gene 1 protein,Zinc finger protein WIG-1,p53-activated gene 608 protein,Pag608,Wig1
Accession No.	Uniprot ID:O54836
Uniprot	O54836
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Zmat3

Note: This product is for in vitro research use only