

# TYRO protein tyrosine kinase-binding protein

Catalog No: #AP89208



Package Size: #AP89208-1 50ug #AP89208-2 100ug #AP89208-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TYRO protein tyrosine kinase-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	DNAX-activation protein 12,DAP12
Accession No.	Uniprot ID:Q9TU45
Uniprot	Q9TU45
GeneID	397405;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TYROBP

Note: This product is for in vitro research use only