

## Tyrosine-protein kinase ZAP-70

Catalog No: #AP89210



Package Size: #AP89210-1 50ug #AP89210-2 100ug #AP89210-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tyrosine-protein kinase ZAP-70
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	70 kDa zeta-chain associated protein,Syk-related tyrosine kinase,Srk,Zap-70
Accession No.	Uniprot ID:P43404
Uniprot	P43404
GeneID	22637;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Zap70

Note: This product is for in vitro research use only