

## Zinc finger homeobox protein 4

Catalog No: #AP89213



Package Size: #AP89213-1 50ug #AP89213-2 100ug #AP89213-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Zinc finger homeobox protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Zinc finger homeodomain protein 4,ZFH-4
Accession No.	Uniprot ID:Q86UP3
Uniprot	Q86UP3
GeneID	79776;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ZFHX4

Note: This product is for in vitro research use only