

E3 ubiquitin-protein ligase TRIM63

Catalog No: #AP89215



Package Size: #AP89215-1 50ug #AP89215-2 100ug #AP89215-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase TRIM63
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Muscle-specific RING finger protein 1,MuRF-1,RING-type E3 ubiquitin transferase TRIM63,Tripartite motif-containing protein 63,Murf1,Rf1
Accession No.	Uniprot ID:Q38HM4
Uniprot	Q38HM4
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Trim63

Note: This product is for in vitro research use only