

## Trem-like transcript 2 protein

Catalog No: #AP89219



Package Size: #AP89219-1 50ug #AP89219-2 100ug #AP89219-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Trem-like transcript 2 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TLT-2, Triggering receptor expressed on myeloid cells-like protein 2, C6orf76, TLT2, UNQ6268/PRO20473
Accession No.	Uniprot ID:Q5T2D2
Uniprot	Q5T2D2
GeneID	79865;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: TREML2

Note: This product is for in vitro research use only