

WAP four-disulfide core domain protein 2

Catalog No: #AP89221



Package Size: #AP89221-1 50ug #AP89221-2 100ug #AP89221-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	WAP four-disulfide core domain protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Epididymal secretory protein E4, Major epididymis-specific protein E4, Putative protease inhibitor WAP5, HE4, WAP5
Accession No.	Uniprot ID:Q14508
Uniprot	Q14508
GeneID	10406;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: WFDC2

Note: This product is for in vitro research use only