

Ubiquitin carboxyl-terminal hydrolase isozyme L1

Catalog No: #AP89223



Package Size: #AP89223-1 50ug #AP89223-2 100ug #AP89223-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ubiquitin carboxyl-terminal hydrolase isozyme L1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	UCH-L1,Neuron cytoplasmic protein 9.5,PGP 9.5,PGP9.5,Ubiquitin thioesterase L1
Accession No.	Uniprot ID:Q00981
Uniprot	Q00981
GeneID	29545;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Uchl1

Note: This product is for in vitro research use only