

Transcription factor HIVEP2

Catalog No: #AP89227



Package Size: #AP89227-1 50ug #AP89227-2 100ug #AP89227-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription factor HIVEP2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Human immunodeficiency virus type I enhancer-binding protein 2 homolog,Myc intron-binding protein 1,MIBP-1,Mibp1
Accession No.	Uniprot ID:Q3UHF7
Uniprot	Q3UHF7
GeneID	15273;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hivep2

Note: This product is for in vitro research use only