

## WD repeat-containing protein WRAP73

Catalog No: #AP89228



Package Size: #AP89228-1 50ug #AP89228-2 100ug #AP89228-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	WD repeat-containing protein WRAP73
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	WD repeat-containing protein 8,WD repeat-containing protein antisense to TP73 gene,Dd57,Wdr8
Accession No.	Uniprot ID:Q9JM98
Uniprot	Q9JM98
GeneID	59002;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Wrap73

Note: This product is for in vitro research use only