

Vascular cell adhesion protein 1

Catalog No: #AP89231



Package Size: #AP89231-1 50ug #AP89231-2 100ug #AP89231-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vascular cell adhesion protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	V-CAM 1,CD106,Vcam-1
Accession No.	Uniprot ID:P29533
Uniprot	P29533
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Vcam1

Note: This product is for in vitro research use only