

von Willebrand factor A domain-containing protein 2

Catalog No: #AP89235



Package Size: #AP89235-1 50ug #AP89235-2 100ug #AP89235-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	von Willebrand factor A domain-containing protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	A domain-containing protein similar to matrilin and collagen,AMACO,Amaco
Accession No.	Uniprot ID:Q70UZ7
Uniprot	Q70UZ7
GeneID	240675;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Vwa2

Note: This product is for in vitro research use only