

Translocator protein

Catalog No: #AP89236



Package Size: #AP89236-1 50ug #AP89236-2 100ug #AP89236-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Translocator protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Mitochondrial benzodiazepine receptor,PKBS,Peripheral-type benzodiazepine receptor,PBR,BZRP,MBR
Accession No.	Uniprot ID:P30536
Uniprot	P30536
GeneID	706;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TSPO

Note: This product is for in vitro research use only