

Uveal autoantigen with coiled-coil domains and ankyrin repeats protein

Catalog No: #AP89237

Package Size: #AP89237-1 50ug #AP89237-2 100ug #AP89237-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Uveal autoantigen with coiled-coil domains and ankyrin repeats protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Beta-actin-binding protein,BetaCAP73
Accession No.	Uniprot ID:Q8HYY4
Uniprot	Q8HYY4
GeneID	281559;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:UACA

Note: This product is for in vitro research use only