

## Ubiquilin-1

Catalog No: #AP89238



Package Size: #AP89238-1 50ug #AP89238-2 100ug #AP89238-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ubiquilin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein linking IAP with cytoskeleton 1, PLIC-1, DA41, PLIC1
Accession No.	Uniprot ID:Q9UMX0
Uniprot	Q9UMX0
GeneID	29979;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: UBQLN1

Note: This product is for in vitro research use only