

E3 ubiquitin-protein ligase UBR1

Catalog No: #AP89239



Package Size: #AP89239-1 50ug #AP89239-2 100ug #AP89239-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase UBR1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	N-recognin-1,RING-type E3 ubiquitin transferase UBR1,Ubiquitin-protein ligase E3-alpha-1,Ubiquitin-protein ligase E3-alpha-I
Accession No.	Uniprot ID:Q8I WV7
Uniprot	Q8I WV7
GeneID	197131;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:UBR1

Note: This product is for in vitro research use only