

## Vascular endothelial growth factor D

Catalog No: #AP89241



Package Size: #AP89241-1 50ug #AP89241-2 100ug #AP89241-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Vascular endothelial growth factor D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	VEGF-D,c-Fos-induced growth factor,FIGF,Figf
Accession No.	Uniprot ID:P97946
Uniprot	P97946
GeneID	14205;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Vegfd

Note: This product is for in vitro research use only