

Zinc transporter 3

Catalog No: #AP89242



Package Size: #AP89242-1 50ug #AP89242-2 100ug #AP89242-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Zinc transporter 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ZnT-3, Solute carrier family 30 member 3, ZNT3
Accession No.	Uniprot ID:Q08E25
Uniprot	Q08E25
GeneID	512803;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SLC30A3

Note: This product is for in vitro research use only