

Ubiquilin-2

Catalog No: #AP89246



Package Size: #AP89246-1 50ug #AP89246-2 100ug #AP89246-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ubiquilin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Chap1,DSK2 homolog,Protein linking IAP with cytoskeleton 2,PLIC-2,Ubiquitin-like product Chap1/Dsk2,N4BP4,PLIC2,HRIHFB2157
Accession No.	Uniprot ID:Q9UHD9
Uniprot	Q9UHD9
GeneID	29978;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:UBQLN2

Note: This product is for in vitro research use only