

Alpha-tocopherol transfer protein

Catalog No: #AP89249



Package Size: #AP89249-1 50ug #AP89249-2 100ug #AP89249-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-tocopherol transfer protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Alpha-TTP, Tpp1
Accession No.	Uniprot ID:P41034
Uniprot	P41034
GeneID	25571;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ttpa

Note: This product is for in vitro research use only