

# X-box-binding protein 1

Catalog No: #AP89251



Package Size: #AP89251-1 50ug #AP89251-2 100ug #AP89251-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	X-box-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	XBP-1
Accession No.	Uniprot ID:Q3SZZ2
Uniprot	Q3SZZ2
GeneID	541236;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: XBP1

Note: This product is for in vitro research use only