

Vasoactive intestinal peptide

Catalog No: #AP89256



Package Size: #AP89256-1 50ug #AP89256-2 100ug #AP89256-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vasoactive intestinal peptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	VIP,Vasoactive intestinal polypeptide
Accession No.	Uniprot ID:P63289
Uniprot	P63289
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:VIP

Note: This product is for in vitro research use only