

Vasoactive intestinal polypeptide receptor 1

Catalog No: #AP89257



Package Size: #AP89257-1 50ug #AP89257-2 100ug #AP89257-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vasoactive intestinal polypeptide receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	VIP-R-1,Pituitary adenylate cyclase-activating polypeptide type II receptor,PACAP type II receptor,VPAC1
Accession No.	Uniprot ID:P32241
Uniprot	P32241
GeneID	7433;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:VIPR1

Note: This product is for in vitro research use only