

Tetraspanin-7

Catalog No: #AP89260



Package Size: #AP89260-1 50ug #AP89260-2 100ug #AP89260-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tetraspanin-7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Tspan-7,Cell surface glycoprotein A15,PE31,TALLA homolog,Transmembrane 4 superfamily member 2,CD231,Mxs1,Tm4sf2
Accession No.	Uniprot ID:Q62283
Uniprot	Q62283
GeneID	21912;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tspan7

Note: This product is for in vitro research use only