

## mRNA decay activator protein ZFP36

Catalog No: #AP89261



Package Size: #AP89261-1 50ug #AP89261-2 100ug #AP89261-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	mRNA decay activator protein ZFP36
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	G0/G1 switch regulatory protein 24, Growth factor-inducible nuclear protein NUP475, Tristetraprolin, Zinc finger protein 36, Zfp-36, G0S24, NUP475, RNF162A, TIS11A, TTP
Accession No.	Uniprot ID:P26651
Uniprot	P26651
GeneID	7538;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: ZFP36

Note: This product is for in vitro research use only