

Thymosin beta-4

Catalog No: #AP89262



Package Size: #AP89262-1 50ug #AP89262-2 100ug #AP89262-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thymosin beta-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	T beta 4,Thyb4,Tmsb4
Accession No.	Uniprot ID:P62329
Uniprot	P62329
GeneID	81814;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tmsb4x

Note: This product is for in vitro research use only