

SUMO-conjugating enzyme UBC9

Catalog No: #AP89263



Package Size: #AP89263-1 50ug #AP89263-2 100ug #AP89263-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | SUMO-conjugating enzyme UBC9 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | RING-type E3 SUMO transferase UBC9,SUMO-protein ligase,Ubiquitin carrier protein 9,Ubiquitin carrier protein I,Ubiquitin-conjugating enzyme E2 I,Ubiquitin-conjugating enzyme UbcE2A,Ubiquitin-protein ligase I,Ubce9 |
| Accession No. | Uniprot ID:P63281 |
| Uniprot | P63281 |
| GeneID | 25573; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ube2i

Note: This product is for in vitro research use only