

## Vasopressin V1b receptor

Catalog No: #AP89265



Package Size: #AP89265-1 50ug #AP89265-2 100ug #AP89265-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Vasopressin V1b receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	V1bR,AVPR V1b,AVPR V3,Antidiuretic hormone receptor 1b,Vasopressin V3 receptor,V1br
Accession No.	Uniprot ID:Q9WU02
Uniprot	Q9WU02
GeneID	26361;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Avpr1b

Note: This product is for in vitro research use only