

Urokinase-type plasminogen activator

Catalog No: #AP89274



Package Size: #AP89274-1 50ug #AP89274-2 100ug #AP89274-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Urokinase-type plasminogen activator
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	U-plasminogen activator
Accession No.	Uniprot ID:P04185
Uniprot	P04185
GeneID	396985;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLAU

Note: This product is for in vitro research use only