

Putative T-complex protein 1 subunit theta-like 1

Catalog No: #AP89277



Package Size: #AP89277-1 50ug #AP89277-2 100ug #AP89277-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Putative T-complex protein 1 subunit theta-like 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CCT8L1
Accession No.	Uniprot ID:A6NM43
Uniprot	A6NM43
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CCT8L1P

Note: This product is for in vitro research use only