

von Willebrand factor

Catalog No: #AP89279



Package Size: #AP89279-1 50ug #AP89279-2 100ug #AP89279-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	von Willebrand factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	vWF,F8VWF
Accession No.	Uniprot ID:P80012
Uniprot	P80012
GeneID	280958;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:VWF

Note: This product is for in vitro research use only