

NEDD8-activating enzyme E1 catalytic subunit

Catalog No: #AP89280



Package Size: #AP89280-1 50ug #AP89280-2 100ug #AP89280-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NEDD8-activating enzyme E1 catalytic subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NEDD8-activating enzyme E1C,Ubiquitin-activating enzyme E1C,Ubiquitin-like modifier-activating enzyme 3,Ubiquitin-activating enzyme 3,Ube1c
Accession No.	Uniprot ID:Q8C878
Uniprot	Q8C878
GeneID	22200;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Uba3

Note: This product is for in vitro research use only