

Zymogen granule membrane protein 16

Catalog No: #AP89289



Package Size: #AP89289-1 50ug #AP89289-2 100ug #AP89289-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Zymogen granule membrane protein 16
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Zymogen granule protein 16,Secretory lectin ZG16
Accession No.	Uniprot ID:Q8CJD3
Uniprot	Q8CJD3
GeneID	171449;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Zg16

Note: This product is for in vitro research use only