

Fractalkine

Catalog No: #AP89291



Package Size: #AP89291-1 50ug #AP89291-2 100ug #AP89291-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fractalkine
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C-X3-C motif chemokine 1,CX3C membrane-anchored chemokine,Neurotactin,Small-inducible cytokine D1,FKN,NTT,SCYD1,A-152E5.2
Accession No.	Uniprot ID:P78423
Uniprot	P78423
GeneID	6376;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CX3CL1

Note: This product is for in vitro research use only