

Vascular cell adhesion protein 1

Catalog No: #AP89293



Package Size: #AP89293-1 50ug #AP89293-2 100ug #AP89293-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vascular cell adhesion protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	V-CAM 1,CD106,Vcam-1
Accession No.	Uniprot ID:P29534
Uniprot	P29534
GeneID	25361;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Vcam1

Note: This product is for in vitro research use only