

Vascular endothelial growth factor receptor 3

Catalog No: #AP89298



Package Size: #AP89298-1 50ug #AP89298-2 100ug #AP89298-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vascular endothelial growth factor receptor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	VEGFR-3,Fms-like tyrosine kinase 4,FLT-4,Tyrosine-protein kinase receptor FLT4,Flt-4,Vegfr3
Accession No.	Uniprot ID:Q91ZT1
Uniprot	Q91ZT1
GeneID	114110;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Flt4

Note: This product is for in vitro research use only