

WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 2

Catalog No: #AP89301

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89301-1 50ug #AP89301-2 100ug #AP89301-3 1mg

Description

Product Name	WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Growth and differentiation factor-associated serum protein 1,GASP-1,WAP, follistatin, immunoglobulin, Kunitz and NTR domain-containing-related protein,WFIKKN-related protein,GASP1,WFIKKNRP,UNQ9235/PRO31996
Accession No.	Uniprot ID:Q8TEU8
Uniprot	Q8TEU8
GeneID	124857;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:WFIKKN2

Note: This product is for in vitro research use only