

Vascular endothelial growth factor receptor 2

Catalog No: #AP89306



Package Size: #AP89306-1 50ug #AP89306-2 100ug #AP89306-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vascular endothelial growth factor receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	VEGFR-2,Fetal liver kinase 1,FLK-1,Kinase NYK,Protein-tyrosine kinase receptor flk-1,CD309,Flk-1,Flk1
Accession No.	Uniprot ID:P35918
Uniprot	P35918
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Kdr

Note: This product is for in vitro research use only