

Upstream stimulatory factor 1

Catalog No: #AP89308

Package Size: #AP89308-1 50ug #AP89308-2 100ug #AP89308-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Upstream stimulatory factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Major late transcription factor 1,Usf
Accession No.	Uniprot ID:Q61069
Uniprot	Q61069
GeneID	22278;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Usf1

Note: This product is for in vitro research use only