

Protein Wnt-4

Catalog No: #AP89312



Package Size: #AP89312-1 50ug #AP89312-2 100ug #AP89312-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein Wnt-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	UNQ426/PRO864
Accession No.	Uniprot ID:P56705
Uniprot	P56705
GeneID	54361;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: WNT4

Note: This product is for in vitro research use only