

Chemokine XC receptor 1

Catalog No: #AP89314



Package Size: #AP89314-1 50ug #AP89314-2 100ug #AP89314-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chemokine XC receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	G-protein coupled receptor 5,Lymphotactin receptor,XC chemokine receptor 1,CCXCR1,GPR5
Accession No.	Uniprot ID:P46094
Uniprot	P46094
GeneID	2829;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:XCR1

Note: This product is for in vitro research use only