

Zinc finger E-box-binding homeobox 1

Catalog No: #AP89315



Package Size: #AP89315-1 50ug #AP89315-2 100ug #AP89315-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Zinc finger E-box-binding homeobox 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Delta EF1,Delta-crystallin enhancer-binding factor,Transcription factor 8,TCF-8,TCF8
Accession No.	Uniprot ID:P36197
Uniprot	P36197
GeneID	396029;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ZEB1

Note: This product is for in vitro research use only