

Zona pellucida sperm-binding protein 3

Catalog No: #AP89317



Package Size: #AP89317-1 50ug #AP89317-2 100ug #AP89317-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Zona pellucida sperm-binding protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Sperm receptor,Zona pellucida glycoprotein 3,Zp-3,Zona pellucida glycoprotein 3-beta,Zp-3-beta,Zona pellucida protein C,ZP3B,ZPC
Accession No.	Uniprot ID:P42098
Uniprot	P42098
GeneID	396927;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ZP3

Note: This product is for in vitro research use only