

Xylosyltransferase 1

Catalog No: #AP89319



Package Size: #AP89319-1 50ug #AP89319-2 100ug #AP89319-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Xylosyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Peptide O-xylosyltransferase 1,Xylosyltransferase I,XT-I,XT1
Accession No.	Uniprot ID:Q86Y38
Uniprot	Q86Y38
GeneID	64131;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:XYLT1

Note: This product is for in vitro research use only