

Zinc finger protein 654

Catalog No: #AP89326



Package Size: #AP89326-1 50ug #AP89326-2 100ug #AP89326-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Zinc finger protein 654
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Melanoma-associated antigen
Accession No.	Uniprot ID:Q8IZM8
Uniprot	Q8IZM8
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ZNF654

Note: This product is for in vitro research use only