

Zinc finger homeobox protein 3

Catalog No: #AP89328



Package Size: #AP89328-1 50ug #AP89328-2 100ug #AP89328-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Zinc finger homeobox protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AT motif-binding factor 1,AT-binding transcription factor 1,Alpha-fetoprotein enhancer-binding protein,Zinc finger homeodomain protein 3,ZFH-3,ATBF1
Accession No.	Uniprot ID:Q15911
Uniprot	Q15911
GeneID	463;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ZFH3

Note: This product is for in vitro research use only