

Mitochondrial uncoupling protein 2

Catalog No: #AP89332



Package Size: #AP89332-1 50ug #AP89332-2 100ug #AP89332-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mitochondrial uncoupling protein 2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | UCP 2,Solute carrier family 25 member 8,UCPH,Slc25a8 |
| Accession No. | Uniprot ID:P70406 |
| Uniprot | P70406 |
| GeneID | 22228; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ucp2

Note: This product is for in vitro research use only