

## Mitochondrial uncoupling protein 3

Catalog No: #AP89333



Package Size: #AP89333-1 50ug #AP89333-2 100ug #AP89333-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mitochondrial uncoupling protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	UCP 3,Solute carrier family 25 member 9,SLC25A9
Accession No.	Uniprot ID:O97649
Uniprot	O97649
GeneID	397116;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:UCP3

Note: This product is for in vitro research use only