

Vang-like protein 2

Catalog No: #AP89334



Package Size: #AP89334-1 50ug #AP89334-2 100ug #AP89334-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vang-like protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Van Gogh-like protein 2
Accession No.	Uniprot ID:P84889
Uniprot	P84889
GeneID	289229;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Vangl2

Note: This product is for in vitro research use only