

V-type proton ATPase catalytic subunit A

Catalog No: #AP89335



Package Size: #AP89335-1 50ug #AP89335-2 100ug #AP89335-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	V-type proton ATPase catalytic subunit A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	V-ATPase 116 kDa isoform a2,Immune suppressor factor J6B7,ISF,Lysosomal H(+)-transporting ATPase V0 subunit a2,ShIF,Vacuolar proton translocating ATPase 116 kDa subunit a isoform 2,Atp6n1b,Tj6
Accession No.	Uniprot ID:P15920
Uniprot	P15920
GeneID	21871;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Atp6v1a

Note: This product is for in vitro research use only